

TRIGO UP
TRIGO

Gebrauchsanleitung
User's Manual

Komponenten:

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Contents

Accessories	2	Safety information	29
Adjusting your TRIGO	3	Intended use.....	29
Checking the seat position	3	STEPS motor option	29
Adjusting the seat position	3	Before your first ride.....	29
Adjusting the seat inclination	3	Observing traffic regulations	29
Adjusting the seat height.....	4	Avoiding muscle strain	29
Adjusting the steering on the TRIGO UP.....	4	Tipping and overturning	29
Adjusting the steering on the TRIGO	5	Braking.....	30
Adjusting the frame length.....	5	Suitable clothing and helmet.....	30
Riding your TRIGO	7	Risk of injury from rotating parts.....	30
Correct tire pressure.....	7	Pedals	30
Shifting	7	Being visible for others	30
Shifting gears with the twist shifter.....	8	Riding at night	30
Braking.....	8	Technical condition	30
Transporting children.....	8	Replacement parts and accessories	30
Parking brake.....	8	Loading the TRIGO.....	30
Parking the TRIGO upright.....	9	Risk of flammability.....	30
Transporting the TRIGO on a car rack.....	9	Description	31
TRIGO with differential	9	What is the TRIGO?.....	31
Transporting the TRIGO	9	Indications for use.....	31
Separate the front boom	10	Contraindications.....	31
Maintenance and repair	11	Adaptability to individual needs.....	31
Maintenance intervals	11	Multiple users	31
Checking wheels for lateral/radial true.....	12	Protection against environmental influences	31
Checking tires.....	12	Cleaning and disinfection.....	31
Changing tires or repairing a flat.....	12	Applicable standards and regulations	31
Checking the disc brakes.....	13	Frame number and bike coding	31
Checking the chain	13	Conditions of warranty and liability	32
Cleaning and oiling the chain.....	13	Inspection pass	33
Checking the chain tube of the TRIGO UP	14	Ergonomical adjustments	35
Cleaning and disinfecting your TRIGO.....	14	Adjustment conflicts	37
Appendix for dealers	15		
Technical data.....	15		
Tightening torques.....	16		
Accepting delivery	16		
Unpacking and checking contents	16		
Assembling the TRIGO.....	16		
Assembling the TRIGO UP	18		
Removing the wheels	19		
Installing the wheels	19		
Replacing the cassette.....	19		
Replacement parts	21		
TRIGO, complete.....	21		
TRIGO UP, complete.....	22		
TRIGO handlebar	23		
Drive shaft bearing.....	24		
Drive assembly, non-driven wheel.....	25		
Drive assembly, TRIGO.....	26		
Drive assembly, TRIGO UP	27		
Assembly parts for Shimano Nexus			
internal gear hub	28		

Zubehör

Accessories

Shimano STEPS Motor
Shimano STEPS motor

Differential
Differential

Faltverdeck
Foldable Fairing

Akku Lichtanlage
Battery lighting kit

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
Visit our homepage for further information

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
Adjusting your TRIGO

The proper adjustment of the TRIGO is important for riding comfort and performance. The TRIGO can be adjusted for use by riders from approx. 1.40 to 2 m (4'7" to 6'7") in height.

Thanks to the adjustable seat, the TRIGO can be adjusted without tools by 20 cm. For larger adjustments an adaptation of the frame and the chain is necessary. Please refer to chapter "Adjusting the frame length", page 5 on how to adjust the frame length.

In addition to the length adjustment, it is also possible to individually adjust the seat height to your own requirements.

To adjust the seat height, please refer to chapter "Adjusting the seat height" on page 4.

Note:

With the separately available extension kit, the TRIGO can also be adapted for persons with a body size over 2,00m (6'7").

Checking the seat position

Sit down on the TRIGO seat and place your feet on the pedals with the balls of your feet on the pedal axles. If one pedal is at the farthest point, the leg should be almost straight and in a comfortable position. If the pedals are too close (e.g. it is not possible to extend the leg) or too far (e.g. the rider cannot easily reach both pedals), then the length of the TRIGO must be adjusted.

Adjusting the seat position

1. Loosen the levers (1) on the front seat mount.
2. Move the seat to the desired position and tighten the previously released levers.
3. Check the seat position. If any changes are necessary, repeat the steps.

Note:

If you want to adjust the length by more than approx. 20 cm, you have to adjust the frame length (see page 5).


Adjusting the seat inclination

1. Loosen the two quick releases (2) on the rear seat attachment.
2. Move the seat to the desired position and tighten the two quick releases (2).
3. Secure the position additionally by clipping in the seat support clips in the highest possible position (3) on both sides.

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

Adjusting the seat height

Adjusting the seat height at the seat

By means of the mounting plate (1) the seat height can be increased by a total of 5cm in two steps. Loosen the front clamping lever and unscrew the rear clamping lever.

Twist the mounting plate to the desired position. Place the square nut in the groove for the hole and screw both clamping levers tight.

Adjusting the seat height at the frame

By adjusting the frame, the seat height can be increased by a further 3cm. To do this, loosen the four bolts (2) on the main frame and move the front boom with the rear screws to the upper position. Make sure that the bolts are in the stop position.

Caution!

It is absolutely necessary to tighten these bolts with a torque of 18 Nm (14 ft-lbs). Riding with loose bolts poses a high risk of accident and damage to the frame.

Adjusting the steering on the TRIGO UP

Check that the handlebar is adjusted so that your arms feel comfortable and your knees do not bump against the handlebar when pedaling. Painful hands, back or neck pain usually indicate an incorrectly adjusted handlebar and this should then be adjusted.

Inclination

Loosen the quick release (3) on the stem and adjust the inclination of the steering rod.

Caution!

Make sure that the lever of the quick release can not touch the front boom at full steering angle (see pictogram).

Length

Loosen the middle quick release (4) and adjust the length of the upper steering rod. Make sure that the upper steering rod is inserted with an insertion depth of at least 10cm.

Alignment

Loosen the quick release clamp (5) on the upper stem and rotate the handlebar to support your hands, allowing your fingers to easily reach the brake and shift levers when the palm is relaxed.

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
Ensuring proper control

If bolts of the steering rod or the steering rod extension are not tightened, you could lose control. Try to move each part and tighten the parts that move or bring the bike to your dealer.

Adjusting the steering on the TRIGO

You can adjust the handlebar width and tilt to your shoulder width and arm length.

1. Unscrew the screws (1) on each side of the handlebar with a 5mm Allen key.
2. Pull out the handlebar ends to widen the handlebar or push the handlebar ends to narrow it.
3. Adjust the inclination of the handlebars.
4. Tighten the screws.

Caution!

The lower ends of the grips must be inserted at least 4 cm into the handlebar. Major risk of accident! If too much pressure is placed on the grips, the center section could break. Make sure that the handlebar grips cannot bump against the wheels or fenders. Major risk of accident! The steering could lock in the handlebar's maximum left or right position.

Adjusting the frame length

If the adjustment range of the seat is extended or the desired seat angle is not reached due to an extreme size adjustment, the frame length can be set individually.


Adjusting the frame length to the rider additionally optimizes the center of gravity of the trike and thereby improves the cornering stability, the turning circle and the riding comfort.

Note:

A frame length change is only possible in combination with a shortening or extension of the chain. In addition to a chain tool, you will also need the appropriate chain links for an extension. The TRIGO comes with an 8-speed chain.

1. Unscrew the four bolts (2) at the rear frame.
2. Open the chain at the master link and remove the chain. With the TRIGO with hub shifting you can push the front boom slightly backwards to take tension off the chain.

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
3. Move the front boom to the desired position. Shorten or lengthen the chain according to the setting. Mount the chain.

**Note:**

Pull the front boom back slightly to close the master link. The chain tension is then adjusted by a slight push of the front boom.

4. Tighten the four Allen bolts (2) on the rear triangle with a 5mm Allen key.

**Caution!**

It is absolutely necessary to tighten these bolts with a torque of 20 Nm (15 ft-lbs). Riding with loose bolts poses a high risk of accident and damage to the frame.

**Note:**

Any restrictions in handlebar movement could lead to accidents with possible injuries.

**Note:**

With the TRIGO, the movements of the handlebar are transferred to the front wheel via a steering transmission rod. After changing the length of the frame, the handlebar must be reoriented.

5. Adjust the steering. To do this, loosen the M6 screw at the clamp (3) of the steering rod. Align front wheel and steering with each other and tighten the M6 screw.

**Caution!**

It is absolutely necessary to tighten these bolt with a torque of 9 Nm (6.6 ft-lbs). Riding with loose bolts poses a high risk of accident and damage to the frame.

6. Adjust the handlebars and the seat (see page 4).
7. Check the position of the brake and shift cables. Check that the handlebar has full freedom of movement by moving it as far as it will go to the left and to the right.

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Riding your TRIGO

Correct tire pressure

The TRIGO should only be ridden with the correct tire pressure. The tire pressure range is specified on the tire sidewalls.

Shifting

Shifting gears allows you to adjust pedaling strength to the terrain and desired speed. The force required per pedal revolution can be decreased by selecting the appropriate gear. In this way, it is much easier to ride up a steep hill in a lower gear than in a high gear. In a lower gear, however, you are slower and have to pedal faster. High gears are best suited for riding at higher speeds on flat terrain or riding downhill. Select a gear that allows you to pedal at more than 70 pedal revolutions per minute. This cadence is optimal for producing sustained power and preventing excessive strain on the knee joints.

In its standard version, the TRIGO is equipped with a 24-speed derailleur system or with a 8-speed Shimano Nexus internally geared hub. With derailleur systems, the gears can only be shifted when pedaling. Hub gears allow you to shift even when the pedals are not in motion.

Note:

Please note the information provided in the enclosed manual from the gear system manufacturer.

Note:

Before riding the TRIGO on public roads, be sure to familiarize yourself with gear shifting in an open area away from traffic. Shifting gears takes practice, which could distract you from navigating road traffic.

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Shifting gears with the twist shifter

Turn the twist shifter (1) on the right handlebar clockwise to shift to a higher gear for faster riding or counterclockwise to shift to a lower gear for slower riding.

Bikes with derailleur systems are also equipped with a twist shifter for the front derailleur. Turn the handle clockwise, then switch to lower gear for slower riding.

Braking

In its standard version, the TRIGO is equipped with mechanical disc brakes.

The brake levers (2) are mounted on the left and right grips of the handlebar. The left wheel brake is controlled with the left brake lever, and the right wheel brake with the right brake lever.

Note:

The disc brakes of the bike are very powerful and have good modulation. However, it is still possible for the wheels to become locked. Before riding the TRIGO on public roads, be sure to familiarize yourself with the braking force by trying a few braking trials in an open area away from traffic.

Transporting children

The special frame of the TRIGO cannot accommodate a child's seat. If you want to transport children, a child trailer must be used. A universal coupling is available as an accessory. For older children who want to help with the pedalling, we recommend the Hase trailer Trets.

Parking brake

The TRIGO is equipped with a parking brake to prevent it from rolling when parked. For this purpose, a push button is provided on a brake lever, with which the brake lever can be fixed. The parking brake is located on the right handlebar handle of the TRIGO UP and on the left handlebar handle of the TRIGO.

1. Press the brake lever and fix it with the push button.
2. The brake is now engaged.
3. Press the lever again. The built-in spring releases the brake.

Note:

Be sure to release the parking brake before you start.

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Parking the trike upright

The TRIGO can be parked upright with the front wheel raised. In this position, it requires less space for parking than a standard bicycle. To bring the TRIGO into the upright parking position, simply lift the front wheel.

Note:

Due to the adjustable seat and the frame length, there is a danger of over-tipping when the settings are extreme. Make sure you have a safe place. Such accidents could cause damage to your TRIGO or other nearby objects, such as bicycles, cars, etc.

Transporting the TRIGO on a car rack

The TRIGO can be transported with a standard roof or rear-mounted bicycle carrier. The rack must be equipped for three bicycles and adjusted so that each TRIGO wheel can be anchored on a separate wheel tray. Secure all wheels and the frame using load straps. Check the mounting to be sure that the TRIGO cannot come loose during the trip. Remove all parts that could fall off in transport.

TRIGO with differential

The optional differential offers even more safety and stability on loose road surfaces. It transfers the power to both rear wheels.

The differential may only be retrofitted by an authorized specialist dealer.

Transporting the TRIGO

The TRIGO can be adjusted to its most compact state for transportation by car, train, or other vehicle (see also page 5). The front boom can be fully inserted and the seat can be removed easily.

1. Before shortening the frame for transport, make a note of the position of the front boom in the main frame. This will make it easier to adjust the TRIGO to the required length after transport.
2. Loosen the quick release clamp of the rear seat attachment and the M8 screw of the front seat attachment. Remove the seat.
3. Tilt the seat rests mounted on the main frame to the front.
4. Slide the seat up to the bend on the front boom.
5. Set the frame to the shortest length. Loosen the steering rod clamp at the TRIGO. You will find a description in chapter "Adjusting the frame length" on page 5.
6. Remove the chain from the front chainring. Be careful not to damage the frame.
7. Loosen the lower quick release (1) on the stem of the

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

TRIGO UP and tilt the handlebar backwards or loosen the clamp on the handlebars of the TRIGO and tilt the handlebars backwards.

8. Loosen the quick release clamp on the fork and remove the front wheel.
9. Now the TRIGO is ready for transport.

**Note:**

Be sure to protect all sensitive parts, such as the disc brakes and rear derailleur, from damage.

The assembly is carried out in reverse order.

Separate the front boom

If the TRIGO is still too long, you can separate the front boom with a few further steps.

Preparation upper handlebar

Open the middle quick release (2) and pull the handlebar out of the lower stem. Wrap the handlebar with padding material and secure it.

Preparation lower handlebar

Loosen the cable on the front derailleur.

Loosen the clamp of the handlebar grips on the handlebar. Disconnect the steering rod with two 10mm wrenches by holding the front steering rod (3) and turning the rear steering rod counter-clockwise. Secure the two sections of the steering rod on the front boom and main frame respectively.

10. Loosen the four M6 Allen bolts (4) on the rear section of the front boom and separate the two frame parts.

If you are planning to take the trike on an airplane or have it shipped by a forwarding company, we recommend packaging it in a cardboard box for protection. Our TRIGO shipping box is very well suited for this purpose.

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Maintenance intervals

Note:

Please have all inspection and maintenance work performed by your authorized HASE retailer.

The first inspection should be done after the initial 300 km (200 miles). All additional inspections should be performed in 3000 km (2000 mile) intervals or at yearly intervals (see also page 33, “Inspection pass”).

Component	Action	Before every ride	Once a month	Once a year
Tires	Check tire pressure Check tread wear and sidewalls	X X		
Brakes	Check pad thickness Check alignment to the brake disc Test braking effectiveness	X	X X	
Brake cables	Check for damage	X		
Chain	Check for wear, oil		X	
Crank bracket	Check bolts			X
Wheels	Check for lateral/radial true Check spoke tension		X X	
Handlebar	Inspect			X
Steering rod	Check for play			X
Headset	Check for play			X
Hubs	Check for play			X
Pedals	Check for axle play			X
Rear derailleur	Clean, lubricate		X	
Derailleur cables	Check for damage		X	
Quick releases	Check for tightness	X		
Bolts and nuts	Check for tightness		X	
Bottom bracket	Check for play			X

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Checking wheels for lateral/radial true

The wheels have been carefully trued prior to shipment. During the break-in period, stress on the wheels may cause the spokes to settle into the rims so that the spoke tension is no longer uniform and the wheel untrue. Have your wheels re-trued by your authorized retailer at the first inspection after approx. 300 km (200 miles). Subsequent tensioning is seldom necessary. Check the spoke tension and lateral/radial true of the wheels regularly.

You can check the spoke tension by squeezing together two adjacent spokes, working all the way around the rim. If you notice major differences in spoke tension, the wheel should be re-trued. You can check whether a rim is true by attaching a wire to the fork or rear frame and positioning it so that it almost touches the rim. Lift the TRIGO and rotate the wheel. If portions of the rim move away from or come in contact with the wire, the wheel must be taken to a qualified bicycle mechanic for re-truing.

Note:

The rear wheels do not have to be removed for the truing process.

Checking tires

Check the tires for damage and tread wear. Tires with cuts, cracks, or worn tread should be replaced. Check tire pressure before each ride and always ride with the pressure specified on the tire sidewalls.

Note:

Never ride with excessive tire pressure (higher than specified on the sidewall). The tire could burst. Risk of accident and injury!

Changing tires or repairing a flat

The front wheel of the TRIGO is mounted with a quick release. Simply release it to remove the front wheel.

You do not have to remove the back wheels for changing a tire or repairing a flat. The tires can be removed and remounted using a standard tire lever with the wheels attached.

When removing and remounting tires, take care not to damage the fenders and check that the rim tape is positioned correctly on the rim.

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Checking the disc brakes

Please note the information provided in the enclosed manual from the brake manufacturer. Have a qualified bicycle mechanic check the pad thickness regularly. If you notice any unusual braking sounds or inadequate braking strength, consult your bicycle mechanic immediately.

Checking the chain

The trike's chain is longer than that of a standard upright bicycle. A long chain is generally subject to less wear. Therefore, it does not have to be changed as often as shorter chains.

Chains on models with hub gears need replacing even less frequently. You can check the condition of the chain using a chain wear gauge from a bike shop.

We recommend that you have a qualified bicycle mechanic check the condition of the chain and, if necessary, replace the chain, and possibly also the chainring and cassette.

Cleaning and oiling the chain

In order to ensure that the chain does not wear out before the end of its service life, it must be cleaned and lubricated on a regular basis. It is also important that the intermediate chains for hub gears be cleaned and lubricated also.

1. Pour a small amount of suitable chain cleaner onto a clean, lint-free cotton cloth, and use it to wipe down the chain. When doing so, slowly turn the crank backwards (i.e. opposite the riding direction). Repeat this procedure with clean patches of the cotton cloth until the chain is clean. Allow the cleaning agent to evaporate for approx. 1 hour before applying a lubricant.

Note:

Never use solvents or pressure cleaners. These could wash the lubricants out of the chain bearings and shorten the life of the chain. Never allow chain lubricant to come in contact with the brake discs or tires.

This could significantly reduce braking performance or, in the case of lubricant on tires, wheel traction.

2. Apply chain oil, grease, or wax evenly to the chain links while slowly rotating the crank.

Note:

Rain can wash off traces of lubricant, releasing it into the environment. We recommend biodegradable lubricants.

3. Rotate a few turns to allow the lubricant to spread.

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4. Wait several minutes to allow the lubricant to penetrate. Then wipe off all excess lubricant using a clean, dry cloth.

**Note:**

If the chain has reached its wear limit, a simultaneous replacement of the cassette is recommended.

Checking the chain guard tube

The purpose of the chain tube (1) is to protect your clothing from chain lubricant. Over time, friction from the moving chain will cause the tube to wear. Check the ends of the chain tube for any signs of wear. If the ends are completely worn, have your bicycle mechanic replace the chain tube.

**Note:**

In order to replace the chain tube, the chain must be opened at the master link.

Cleaning and disinfecting your TRIGO

Be sure to clean your TRIGO on a regular basis. Dirt and grime, as well as salt from winter roads, can cause damage. To clean your TRIGO, use a sponge and soapy water (ideally with a biodegradable household cleaning product). Stubborn marks on metal surfaces with a matt finish can be removed best by using a sponge with a wheel-rim cleaner from an automotive-supply shop.



You can also clean your TRIGO with plain water. If you use a hose, be sure to spray the water gently. Never use pressure washers. The high pressure could force water into the bearings, damaging them. Wipe all surfaces with a moist, lint-free cloth.

Apply a suitable spray wax or a similar corrosion inhibitor to all surfaces of the TRIGO except for the brake pads, brake discs, grips, brake/gear levers, chain, seat cover, and tires. Then clean and lubricate the chain. The handlebar grips and seat can be cleaned and/or disinfected using standard cleaning agents and disinfectants. Be sure to use a disinfectant that is suitable for plastic materials.

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Appendix for dealers

Technical data

Total length: 179-210 cm (70½"-83")
Minimum length without front wheel: 155 cm (61")
Total width: 83 cm (33")
Total height: 98-106 cm (38½"-42") (Without seat: 56 cm / 22")
Bottom bracket height: 33 cm (13")
Seat height: low seat position 58 cm (23") high seat position 61-66 cm (24"-26")
Weight TRIGO: 20.9 kg (46.1 lbs) Weight TRIGO UP: 22.7 kg (50.0 lbs)
Max. load: 140 kg (308 lbs)
Brake system: mech. disc brakes
Gears TRIGO: 24-speed Shimano derailleur system / 8-speed Shimano Nexus Gears TRIGO UP: 8-speed Shimano Nexus
Main frame: aluminum, powder-coated Front boom: aluminum black anodized
Wheels: 20"
Rims: aluminum
Tires: 58-406 (20 x 2.35")
Fits riders from 1.40 to 2.00 m (4'7" to 6'7") in height

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
Tightening torques

Attention!

Use a torque wrench for all screws. Please check the safety-relevant screws and observe the following table with torques.

Bolt	Tightening torque
M8 lens head bolts at the front boom	13.3 ± 0.7 ft-lbs
M8 Allen bolt in the rear hub	14.7 ± 0.7 ft-lbs
M8 lens head bolts at the fork	13.3 ± 0.7 ft-lbs
M6 outer sunk bolt at the frame separation	5.9 ± 0.7 ft-lbs
M8 counter nut at the steering rod	11 ± 0.7 ft-lbs

The following table lists the standard values for Allen bolts of 8.8 quality (Grade 5)

Bolt	Tightening torque
M 5	4 ft-lbs
M 6	6.6 ft-lbs
M 8	15.5 ft-lbs

Accepting delivery

If you notice the packaging has been damaged, please notify the shipping company directly.

Unpacking and checking contents

Check the contents of the package with your order. In addition to the front boom, rear frame, and the seat, the package should also include: Pedals, chain, front-wheel quick-release, four M6x12 Allen bolts, TRIGO User's Manual as well as manuals from other manufacturers for various components such as gear system and brakes. In addition, the chain guard tubes are included in the TRIGO UP.

Note:

A detailed assembly video for the TRIGO can be found at www.hasebikes.com.

Assembling the TRIGO

1. Remove the top box from the package.
2. Remove the front wheel and the seat from the package.

Note:

Leave the main frame in the package and use it as a mounting aid for the front boom.

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3. Insert the front section of the front boom (1) into the rear section of the front boom (see top photo) and assemble the front boom with the four enclosed M6x12 screws. Use medium strength loctite.
4. Install the front wheel and the pedals.
5. Tilt the TRIGO together with the package to the ground and pull the TRIGO out of the package (see bottom photo).
6. Mount the chain.
7. Slide the front boom to the desired position and tighten the four M8 screws on the main frame.

**Note:**

Note that the enclosed chain is designed for an average body size. Due to the slidable seat, the TRIGO can be adjusted in this setting to a body size of approx. 1,65-1,85 m (5'5"-6'1"). Shorten or lengthen the chain for a different setting.

8. Loosen the rear steering rod at the steering rod clamp (2). Tighten the steering rod with two 10mm wrenches by holding the front steering rod (3) and turning the rear steering rod clockwise.
9. Align the steering and attach the steering rod to the steering rod clamp (2).
10. Mount the seat. Slide the rear seat attachment over the rear seat tubes (4). Loosen the M8 screw (5) on the front seat attachment and tighten the seat to the seat plate using the M8 screw.
11. Adjust the position and angle of the seat. Secure the position additionally with the seat support clips (see page 4).
12. Adjust the handlebar width and tilt (see page 5).
13. Install front derailleur at the front derailleur support. Adjust the gear system. Also observe the enclosed manual of the gear system manufacturer.
14. Lay the cable under the seat with plenty of slack, and attach the gear cable to the front boom in the respective cable guides. Check the cables when the front wheel is turned all the way to one side.
15. Correct the air pressure in the tires.
16. Check all screws for tight fit.
17. Mount any accessories that you have ordered.

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**Note:**

A detailed assembly video for the TRIGO UP can be found at www.hasebikes.com.

Assembling the TRIGO UP

1. Remove the top box from the package.
2. Remove the front wheel and the seat from the package.

**Note:**

Leave the main frame in the package and use it as a mounting aid for the front boom.

3. Insert the front section of the front boom (1) into the rear section of the front boom (see top photo) and assemble the front boom with the four enclosed M6x12 screws.
4. Install the front wheel and the pedals.
5. Tilt the TRIGO UP together with the package to the ground and pull the TRIGO UP out of the package (see bottom photo).
6. Mount the handlebar to the stem.
7. Slide the front boom to the desired position.
8. Mount the chain with the chain guards.
9. Tension the chain by slightly pushing the front boom and tighten the four M8 screws on the main frame

**Note:**

Note that the enclosed chain is designed for an average body size. Due to the slidable seat, the TRIGO UP can be adjusted in this setting to a body size of approx. 1,65-1,85 m (5'5"-6'1").

Shorten or lengthen the chain for a different setting.

10. Mount the seat. Slide the seat attachment over the rear seat tubes (4). Loosen the M8 screw (5) on the front seat attachment and tighten the seat to the seat plate using the M8 screw.
11. Adjust the position and angle of the seat. Secure the position additionally with the seat support clips (see page 4).
12. Adjust the handlebar (see page 5).
13. Correct the air pressure in the tires.
14. Fix the brake and gear cables in the cable guides provided for this purpose. Check the cables when the front wheel is turned all the way to one side.
15. Check all bolts for tight fit.
16. Mount any accessories that you have ordered.

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
**Removing the wheels**

1. Loosen the brake (1) on the adapter.
2. Hold the wheel firmly and turn the M8 screw (2) counter-clockwise using a 6 mm Allen key.
3. Remove the wheels together with the brake disc.

Installing the wheels

1. Turn the axle so that the feather key (3) points upwards. Position the groove (4) in the hub to the feather key and slide the wheel together with the brake disc onto the axle.
2. Tighten the wheel with the washer and the M8 screw (1).

**Attention!**

It is absolutely necessary to tighten these bolts with a torque of 21 Nm (17 ft-lbs). Use medium strength loctite. Riding with loose bolts poses a high risk of accident and damage to the frame.

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

**Replacing the cassette**

When replacing the chain, check the cassette for wear. Often, a new chain can not be fitted with a worn cassette. To replace the cassette, the right axle must be pushed out to the side.

Removing the cassette

1. Remove the chain from the derailleur.
2. Loosen the grub screws at the adjusting collar (1) on the right axle.
3. Relax the flange bearings by loosening the four self-locking M6 nuts (2).
4. Hold the right wheel and loosen the M8 screw (3) with a 6mm Allen key from the wheel hub.
5. Slide the right axle out so that you can pull out the cassette along with the freewheel from the axle. Use a slide hammer if the axle can not be pulled out by hand.

**Note:**

To avoid damage to the drive shaft, the slide hammer thread rod must be screwed into the thread as far as it will go. Alternatively, the brake caliper can be removed and the wheel used as a pulling aid.

6. Now the cassette can be dismantled and replaced.

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
**Installing the drive shaft**

7. Lubricate the sliding surfaces in the cassette carrier (5) and the brass bushing (6) in the freewheel body generously with grease.
8. Place the freewheel with the new cassette on the cassette support (5).
9. Push the drive shaft through the brass bushing (6) into its end position. Observe the orientation of the feather key (7) to the hub driver and the joint on the shaft to the drive attachment (4).
10. Tighten the M8 screw (3) in the wheel hub.

**Attention!**

It is absolutely necessary to tighten these bolts with a torque of 21 Nm (15.5 ft-lbs). Use medium strength loctite. Riding with loose bolts poses a high risk of accident and damage to the frame.

11. Tighten the four M6 self-locking nuts (2) to the flange bearings.
12. Press the adjusting collar (1) firmly against the outer flange bearing and tighten the two grub screws.
13. Mount the chain and adjust the gear system and the brakes if necessary.

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
Position	Description	Art. No.	Number
1	Seat cover	25744	1
2	Handlebar grips	-	1
3	Bottom bracket	25742	1
4	Front derailleur support	25793	1
5	Steering rod support	25792	1
6	Steering rod front/rear	25757/25756	1
7	20" front wheel	24180	1
8	20" fork	26004	1
9	Seat stay 18mm	25747	1
10	Parking brake, mechanical	-	1
11	Seat stay 22mm	25748	2
12	Seat stay diagonal	25749	2
13	20" rear wheel	24264	1
14	Seat support plate	25746	1

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
Position	Description	Art. No.	Number
1	Seat cover	25744	1
2	Handlebar	-	1
3	Upper stem	26125	1
4	Lower stem	25791	1
5	20" front wheel	24180	1
6	20" fork	26005	1
7	Chain guard tube	-	2
8	Seat stay 18mm	25747	2
9	Parking brake, mechanical	-	1
10	Seat stay 22mm	25748	2
11	Seat stay diagonal	25749	2
12	20" rear wheel	24264	2
13	Bottom bracket	25742	1
14	Seat support plate	25746	1

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

Position	Description	Art. No.	Number
1	M6x45 cup head square neck bolt	-	2
2	M6x25 cup head square neck bolt	-	6
3	Flanged bearing Ø17mm	25765	3
4	Bearing distance shell	25751	2
5	Seat stay 18mm	25747	2
6	M6 self-locking nut	20643	8
7	M6 washer D=2xd	20581	6
8	Spherical bearing	-	1
	Bearing housing	25763	2
	Spherical bearing shell	25762	1
	Grooved ball bearing 6805	25764	2

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

Position	Description	Art. No.	Number
1	Attachment with hexagon	25770	1
2	Cassette support	25787	1
3	8 speed cassette with freewheel body	-	1
4	Brass bushing for freewheel	25785	1
5	Allen bolt M5x10	20687	7
6	Derailleur hanger	25788	1
7	Ball bearing with grub screws, cylindrical	25781	1
8	Drive shaft TRIGO	25759	1
9	Feather key, VA 5x5x30, DIN 6885	26124	1
10	Adjusting collar with grub screws	25760	1
11	Brake disc 160mm	21595	1
12	TRIGO hub	26002	1
13	M8 washer, D=3d	20628	1
14	Allen bolt M8x12	20813	1

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Safety information

Intended use

The TRIGO is a recumbent trike intended for use on paved roads and bicycle paths, as well as on roads or trails surfaced with sand, gravel, or similar materials (e.g. forest paths, dirt roads). Your safety on roads, trails, and off-road terrain will depend on your riding speed. The higher your speed, the greater your risk! The use of the TRIGO for competition or any form of racing is strictly forbidden.

Period of use

With proper use, taking into account the instructions for use and depending on the frequency of use, the expected useful life is up to eight years.

STEPS motor option

The frame of the TRIGO has been tested in combination with the optionally available Shimano STEPS pedal-assist system and satisfies the requirements of the safety standard EN 15194 for EPACs. All Hase Bikes vehicles purchased with a Shimano STEPS pedal-assist system come with user manuals from the manufacturer Shimano, which include all the information you need for the safe use of the Shimano STEPS pedal-assist system. Upon receipt of your Hase Bikes vehicle, please check to make sure that the current Shimano user manuals on the Shimano STEPS E6000 series and the Shimano eTube project are included in the scope of delivery.

Before your first ride...

Your authorized retailer has assembled your new TRIGO and adjusted it to your size requirements. Please refer to the section “Adjusting your TRIGO” on page 4 if you need to make adjustments for another rider. Before setting out on your first trip, please take some time to familiarize yourself with the TRIGO. Due to its unique construction, there are significant differences in handling in comparison with a standard bicycle. Before riding the TRIGO on public roads, it is important to familiarize yourself with the handling in a quiet, open area, away from traffic.

The maintenance measures required before each ride are to be performed by the user (individuals with no physical, mental, or visual impairments).

Observing traffic regulations

The TRIGO is a piece of sports equipment and is delivered in its standard form without a lighting system. Therefore, it cannot be used on public roads until the required accessories have been installed. Consult your authorized retailer for information on the accessories required. All equipment should be installed by a qualified bike mechanic. If you plan on installing or making changes to the lighting system yourself, please study the applicable regulations before you begin. Traffic regulations are different in every country. It is important to choose a lighting system, reflectors, and a bell in accordance with the applicable regulations of the country in which the TRIGO will be used. We offer accessories that comply with the German road traffic regulations (StVZO).

Avoiding muscle strain

In a recumbent position, the rider uses different muscle groups than with a standard bicycle. Riders with little or no experience on recumbent bikes/trikes should begin with short, easy trips in order to avoid excessive muscle strain.

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Tipping and overturning

Recumbent tricycles have a low center of gravity. This allows for cornering at higher speeds. The centrifugal force felt when turning at high speeds, however, can cause the TRIGO to skid or even tip over, which could cause injury to the rider or others. Avoid making hectic steering movements or braking in curves. This increases the TRIGO's tendency to tip. Be careful when stopping on inclines. The TRIGO could roll backwards and, if brakes are applied abruptly, even overturn. Always wear a helmet for your own safety.

Braking

The disc brakes of the TRIGO are very powerful and have good modulation. However, it is still possible for the wheels to become locked. Before your first trip in public traffic, be sure to take the time to familiarize yourself with the braking strength by trying a few braking trials in an open area away from road vehicles.

Suitable clothing and helmet

Do not wear loose clothing that could get caught in the wheels or gears. Beware of scarves and clothing strings/cords that hang down. Always wear a helmet in order to avoid head injury in the case of an accident.

Risk of injury from rotating parts

If a patient seated in a TRIGO that is being pulled as a trailer allows his/her hands or fingers to hang down beside or below the seat, there is a risk that they could get caught in the disc brakes and/or wheels. To prevent the associated injuries, it is essential to ensure that patients keep their hands/fingers on the handlebar grips or above the seat at all times!

Pedals

In a recumbent position, it is especially easy for the feet to slip from the pedals. Always wear shoes with good tread. Maintain sufficient pressure on the pedals. Use toe clips and straps (slip protection) or clipless pedals.

Being visible for others

The TRIGO is lower than an upright bicycle. This may make it more difficult for other road vehicles to see you. Always ride with a flag. It attracts the attention of other road users.

Riding at night

Always use lights when riding at night or in conditions of low visibility, such as dusk and dawn.

Technical condition

The TRIGO should only be ridden in perfect technical condition. Check for damage before each trip. Take the TRIGO to your bike shop for regular maintenance. The recommended maintenance intervals can be found in the section "Maintenance intervals" on page 11.

Replacement parts and accessories

Use only the replacement parts and accessories recommended by HASE BIKES for the TRIGO. If other parts are used, the warranty will lose its validity. HASE BIKES is not liable for any resulting damage. Consult your bike shop for more information.

Loading the TRIGO

Be sure not to exceed the maximum permissible combined weight for rider and luggage: 140 kg (308 pounds). All luggage should be secured using the mount behind the seat or a lowrider rack (accessory). Always ensure that the load is distributed evenly.

Risk of flammability

The upholstery material of the seat cover is not resistant to ignition sources as specified in EN 1021-2. Therefore, be sure to keep the seat away from open flames, heat sources, etc.

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Description

What is the TRIGO?

The TRIGO is a recumbent trike and an approved medical device as per the European Directive (2017/745) MDR (Medical Device Directive). The TRIGO is driven by leg power via the pedals, which are specially designed to prevent the feet from slipping out when riding. It is available with an 8-speed internal-gear hub or derailleur gears. On models with the optional differential, both rear wheels are driven, and on standard models, only one rear wheel is driven. The padded seat has been specially designed and constructed with flexible materials to cushion the ride and offer maximum comfort. The TRIGO has a reliable brake system with two independently operated disc brakes.

Indications for use

The TRIGO is suitable for riders with neuromuscular disorders, such as cerebral palsy or muscular dystrophy, who – owing to special physical conditions – are unable to ride standard bicycles or vehicles, even with training wheels. Other indications include coordination disorders, balance disorders, and restricted movement associated with spasticity, paraplegia, dysmelia, amputations, etc.

Contraindications

The TRIGO is not suitable for riders with visual impairment.

Adaptability to individual needs

The TRIGO can be equipped to accommodate riders with various physical impairments. The following are just a few examples of available accessories: single-hand controls, crank shorteners, special pedals with or without calf support, pedal pendulums, lap and shoulder belts, and holders for walking sticks.

Multiple users

TRIGOs that have been used in accordance with their intended use can be reused by a new rider after the TRIGO has undergone an inspection by an authorized retailer.

Protection against environmental influences

The TRIGO is made from special materials that protect it against corrosion: powder-coated aluminum frame parts, plastics, and stainless steel.

Cleaning and disinfection

The handlebar grips and seat can be cleaned and disinfected using household cleaning agents and disinfectants. For more information, see the section “Cleaning and disinfecting your TRIGO” on page 14.

Applicable standards and regulations

The TRIGO has been designed and manufactured in compliance with the following regulations and standards:

European Directive (2017/745) MDR (Medical Device Directive)

Frame number and bike coding

The frame number can be found on the end of the main frame (see illustration on the inside cover of this manual). This is also the only place where bike coding can be engraved. Bike coding helps authorities locate the vehicle owner in the case of theft.

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Conditions of warranty and liability

The statutory warranty applies. The dealer must sell you the HASE BIKE in a safe and roadworthy condition. In addition, the company HASE BIKES offers a warranty of three years (from the date of purchase) against frame breakage in accordance with the following conditions. By registering your bike at <https://hasebikes.com/209-1-Bicycle-passport.html> you can extend this warranty to 5 years.

- HASE BIKES will only be held liable for defects resulting from faulty frame manufacturing and assembly. For assembly parts, the relevant manufacturer is liable.
- This additional warranty only applies to the original owner and to usage in accordance with the user's manual.
- Within the warranty period of 3 years, HASE BIKES will repair, or if necessary replace, any broken frame part free of charge. Any defective parts replaced by HASE BIKES become the property of HASE BIKES.
- Any further claims beyond the scope of this warranty are excluded. HASE BIKES will not cover the costs of assembly or disassembly (e.g., by the retailer) or shipping. We will not accept postal deliveries with insufficient postage.
- The maximum permissible load (rider and luggage) of 140kg (308 lbs) specified for the TRIGO must not be exceeded.
- If the TRIGO is used in competition or subjected to abnormally excessive stress, the statutory and additional warranty becomes void.
- Damage caused by accidents is not covered by the warranty.
- The warranty becomes void if the TRIGO has been fitted with incompatible accessories or if accessories are installed incorrectly. Therefore, always have components replaced and installed by an authorized retailer.
- The warranty becomes void if improper modifications or alterations (e.g., grinding, drilling, bending, etc.) are made to the frame.
- It is not possible to extend the warranty or commence a new warranty period.
- Warranty claims can only be handled by the authorized HASE BIKES retailer from whom the TRIGO was purchased or by one of our distribution partners.
- The TRIGO may only be ridden by children under parental supervision.

Further information

If you have any problems or questions about your therapeutic bike or trike, please contact your authorized dealer or HASE BIKES for assistance.

Any serious incident* that occurs in connection with the product must be reported to the manufacturer and the responsible competent authorities of the member state in which the user and/or patient resides.

*A "serious incident" is any incident that had, could have had, or could have one of the following consequences, either directly or indirectly: a) the death of a patient, user, or other person, b) the serious temporary or permanent deterioration of the health condition of a patient, user, or other person, c) a serious risk to public health.

Disposal instructions

At the end of its useful life, the product should be disposed of properly, in accordance with the local waste disposal regulations (i.e. not as household waste).

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Have your vehicle inspected by your bicycle dealer at regular intervals. During the high season (March through September), it is best to schedule any required repairs or inspections in advance in order to avoid unnecessarily long downtimes and waiting periods.

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Checklist:

- Check pedals for tightness
- Check headset
- Check handlebar for tightness
- Check front and rear brakes
- Check front and rear wheel bearings
- Check bottom bracket
- Check cranks for tightness
- Check seat for tightness
- Check chain tension and wear
- Check tires and air pressure
- Check rims and spokes
- Check all nuts and bolts for tightness
- Check gear system
- Check accessories like motor system

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After 4000 km (2,500 miles) or two years after the date of sale.

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After 8000 km (5,000 miles) or four years after the date of sale.

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After 6000 km (3,700 miles) or three years after the date of sale.

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After 10.000 km (6,200 miles) or five years after the date of sale.

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