

RACER® STROLLERS FAMILY

Special needs strollers



RACER®



RACER EVO®



RACER HOME®

**This is a medical device.
Use it in accordance with the user manual or the label.**

Description

RACER®

Special needs stroller **RACER®** is used for children whose movement possibilities are limited. It may be used in the rehabilitation centre and at user's home. Stable construction perfectly stabilizes the user and it gives him a sense of security. What is more, it enables the correct sitting position and the child has a chance to rest in the reclining position.

Parent is able to take the child to rehabilitation and go for a walk thanks to which the child has the possibility of spending time outdoor. There is no doubt that this special stroller makes life easier for parent regarding care and feeding and, moreover, the child can actively spend time with his family. It is equipped with a vast array of adjustment thanks to which the user feels comfortable in the device. Additionally, the transport of this buggy is easy due to collapsible construction. Concerning upholstery, it is made of breathable and washable fabrics.

Indications

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.





Special features



Stable structure

The device ensures adequate stability and safety of user.



Easy to use

Intuitive way of adjustment and easy assembly.



Reflective elements

The product has reflective elements to ensure your safety.

24⁺

24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.



Colors

ATMO™ upholstery



HYDRO™ upholstery



*Available while stocks last.



Two-piece upholstery

Two-piece upholstery – much more comfortable and durable. Simple disassembly makes it easy to clean and maintain hygiene.



Side supports

Stabilize the position of the user's trunk in the stroller; the upholstery side pocket system allows for their easy removal.



Footrest


Footrest with 2 points feet stabilizing belts.



Standard equipment




Important measurements



Sizes	Seat width	Seat depth	Backrest width	Headrest height	Backrest height	Seat tilt angle*	Footer length	Indicative user height**	Max user weight
Size 1	35 cm	34 cm	18 cm (with side supports) 30 cm (without side supports)	58 cm	58 cm	15°	23 cm	< 120 cm	35 kg
Size 2	36 cm	42 cm	20 cm (with side supports) 32 cm (without side supports)	58 - 62 cm	62 cm	15°	23 cm	< 150 cm	50 kg
Size 3	39 cm	42 cm	20 cm (with side supports) 36 cm (without side supports)	72 - 80 cm	80 cm	15°	25 cm	< 180 cm	75 kg
Size 4	46 cm	45 cm	25 cm (with side supports) 41 cm (without side supports)	72 - 80 cm	80 cm	15°	25 cm	< 180 cm	90 kg

Weight and dimensions

Size	Width	Width after folding	Length	Length after folding	Height	Height after folding	Weight***
Size 1	62 cm	62 cm	102 - 125 cm	103 cm	103 cm	43 cm	24 kg
Size 2	64 cm	64 cm	118 - 140 cm	110 cm	105 cm	54 cm	27 kg
Size 3	70 cm	70 cm	135 - 162 cm	125 cm	125 cm	58 cm	29,4 kg
Size 4	76 cm	76 cm	135 - 162 cm	125 cm	125 cm	58 cm	30 kg



Scan the QR code to see the full table.

*Value measured with the seat loaded.

**Height cannot be the main parameter while choosing the product size.

***The given weight includes a frame with wheels and upholstery. The weight may vary depending on the additional accessories the stroller can be equipped with.

Note: Measurement tolerance range is +/- 1 cm / 1/2 in.
Note: Approximate values, depending on the product configuration.

Description

RACER EVO®

RACER EVO® stroller is equipped with swivel front wheels, thanks to which parents do not have a problem with using this device. This change is really beneficial due to the fact the buggy does not have to be lifted in order to overcome uneven grounds.

Moreover, it is really good for the child because he feels comfortable and, more to the point, it does not disturb the child's perception. Young user has bigger chance to regain fitness and acquire new skills when nothing disturbs the child's perception in the stroller.

Furthermore, front wheels are equipped with fixed direction of movement which is handy regarding uneven surfaces. New footrest adjustment facilitates the use of this rehab stroller and modern design does not suggest that it is created for people who need special care.

Indications

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.





Special features



Stable structure

The device ensures adequate stability and safety of user.



Easy to use

Intuitive way of adjustment and easy assembly.



Reflective elements

The product has reflective elements to ensure your safety.

24⁺

24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.



Colors

ATMO™ upholstery



HYDRO™ upholstery



*Available while stocks last.



Two-piece upholstery

Two-piece upholstery – much more comfortable and durable. Simple disassembly makes it easy to clean and maintain hygiene.



Side supports

Stabilize the position of the user's trunk in the stroller; the upholstery side pocket system allows for their easy removal.



Footrest


Footrest with 2 points foot stabilizing belts.



Standard equipment




Important measurements



Sizes	Seat width	Seat depth	Backrest width	Headrest height	Backrest height	Seat tilt angle*	Footer length	Indicative user height**	Max user weight
Size 1	35 cm	34 cm	18 cm (with side supports) 30 cm (without side supports)	58 cm	58 cm	15°	23 cm	< 120 cm	35 kg
Size 2	36 cm	42 cm	20 cm (with side supports) 32 cm (without side supports)	58 - 62 cm	62 cm	15°	23 cm	< 150 cm	50 kg
Size 3	39 cm	42 cm	20 cm (with side supports) 36 cm (without side supports)	72 - 80 cm	80 cm	15°	25 cm	< 180 cm	75 kg
Size 4	46 cm	45 cm	25 cm (with side supports) 41 cm (without side supports)	72 - 80 cm	80 cm	15°	25 cm	< 180 cm	90 kg

Weight and dimensions

Size	Width	Width after folding	Length	Length after folding	Height	Height after folding	Weight***
Size 1	64 cm	64 cm	105 - 125 cm	102 cm	101 cm	48 cm	22 kg
Size 2	65 cm	65 cm	123 - 146 cm	110 cm	104 cm	52 cm	24,8 kg
Size 3	67 cm	67 cm	135 - 162 cm	128 cm	122 cm	45 cm	28 kg
Size 4	74 cm	74 cm	135 - 162 cm	128 cm	118 cm	45 cm	28,6 kg



Scan the QR code to see the full table.

*Value measured with the seat loaded.

**Height cannot be the main parameter while choosing the product size.

***The given weight includes a frame with wheels and upholstery. The weight may vary depending on the additional accessories the stroller can be equipped with.

Note: Measurement tolerance range is +/- 1 cm / 1/2 in.
Note: Approximate values, depending on the product configuration.

Description

RACER HOME®

RACER HOME® another product from the family of the special needs strollers **RACER®**. It is the response for our customers' individual needs. It is an indoor version of a **RACER®** stroller that is well known and appreciated by our users and their caretakers.

RACER HOME® allows transporting children and adults with limited locomotor abilities. It is perfect both for private apartments and rehabilitation centers. The sizes 1 – 3 easily fit into 60 cm wide door, which is especially useful in case of narrow rooms and elevators.

A stable construction and a system of belts and supports make the user feel safe and secure.

The stroller retains all the advantages and stabilization systems of the standard **RACER®** stroller as well as the possibility of mounting all the existing accessories.

RACER HOME® has been made of durable materials which will allow you to use it for many years. The stroller can be folded and unfolded very quickly and easily, that is why transporting it in a car or storing it does not cause any problems.

Indications

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.





Special features



Stable structure

The device ensures adequate stability and safety of user.



Easy to use

Intuitive way of adjustment and easy assembly.



Reflective elements

The product has reflective elements to ensure your safety.

24⁺

24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.



Colors

ATMO™ upholstery



Grass
Flex upholstery BLACK



Sea
Flex upholstery BLACK

HYDRO™ upholstery



Lime
Flex upholstery STEEL*



Aqua
Flex upholstery STEEL*



Steel
Flex upholstery STEEL*



Cobalt
Flex upholstery STEEL*



Denim*
Flex upholstery STEEL*

*Available while stocks last.



Two-piece upholstery

Two-piece upholstery – much more comfortable and durable. Simple disassembly makes it easy to clean and maintain hygiene.



Side supports

Stabilize the position of the user's trunk in the stroller; the upholstery side pocket system allows for their easy removal.



Quiet swivel wheels


Four very quiet swivel wheels, each with a brake, make it easier to maneuver the stroller in narrow places. Additionally, the wheels are made of light rubber which does not leave any marks on the floor.



Standard equipment



Important measurements



Sizes	Seat width	Seat depth	Backrest width	Headrest height	Backrest height	Seat tilt angle*	Footer length	Indicative user height**	Max user weight
Size 1	35 cm	34 cm	18 cm (with side supports) 30 cm (without side supports)	58 cm	58 cm	15°	23 cm	< 120 cm	35 kg
Size 2	36 cm	42 cm	20 cm (with side supports) 32 cm (without side supports)	58 - 62 cm	62 cm	15°	23 cm	< 150 cm	50 kg
Size 3	39 cm	42 cm	20 cm (with side supports) 36 cm (without side supports)	72 - 80 cm	80 cm	15°	25 cm	< 180 cm	75 kg
Size 4	46 cm	45 cm	25 cm (with side supports) 41 cm (without side supports)	72 - 80 cm	80 cm	15°	25 cm	< 180 cm	90 kg

Weight and dimensions

Size	Width	Width after folding	Length	Length after folding	Height	Height after folding	Weight***
Size 1	50 cm	50 cm	101 - 122 cm	90 cm	101 cm	47 cm	20 kg
Size 2	53 cm	53 cm	120 - 143 cm	100 cm	99 cm	46 cm	22,6 kg
Size 3	56 cm	56 cm	130 - 156 cm	130 cm	120 cm	45 cm	25,8 kg
Size 4	62 cm	62 cm	132 - 158 cm	130 cm	122 cm	45 cm	26,8 kg



Scan the QR code to see the full table.

*Value measured with the seat loaded.

**Height cannot be the main parameter while choosing the product size.

***The given weight includes a frame with wheels and upholstery. The weight may vary depending on the additional accessories the stroller can be equipped with.

Note: Measurement tolerance range is +/- 1 cm / 1/2 in.
Note: Approximate values, depending on the product configuration.

Accessories RACER® / RACER EVO® / RACER HOME®



Narrowing pads 6/10 cm

Made of elastic foam trimmed with soft material. They are attached to the seat upholstery from both sides.



Pelvis stabilizer

Prevents the pelvis and hips from sliding forward, protecting the user from slipping out of the seat. It helps maintain proper sitting position, which is especially important for people with increased muscle tension.



Elastico cushion for seat

It precisely adapts to the patient's body. Made of memory foam covered with breathable and spatial fabric. It is indicated in anti-bedsore prophylaxis.




Folding canopy

Protects against both sun and rain.



Mosquito net

Excellent protection against insects



Scan the QR code, to see all the accessories.