Product Code	APH041
Bed Dimensions	235 x 105 cm
Mattress Platform Dimensions	200 x 89 cm
Mattress Dimensions	Backrest section: 89 x 72 x 13 cm Thigh rest section: 89 x 66 x 13 cm Calf rest section: 89 x 31 x 13 cm Footrest section: 89 x 20 x 13 cm
Bed Weight (excluding side rails and mattress)	108 kg
Mattress Weight	11 kg
Maximum Weight Limit	28.3 stone / 180 kg
Risk Category	High Risk
Maximum Height Position	72.5 cm
Lowest Height Position	40.5 cm
Height Adjustment	32 cm
Headboard Height from mattress platform	48 cm
Angle of Backrest	72°
Angle of Thigh rest	25°
Angle of Calf rest (after Thigh rest is raised)	3°
Angle of Anti-Trendelenburg tilt	11°
Motor Brand	TiMOTION
Electrically operated backrest, knee break, height, rotate and tilt	✓
Auto-regression function	✓
Innovative lateral support backrest	✓

Due to ongoing research and development, Apollo Healthcare Technologies Ltd, reserve the right to change specifications without prior notice.

If you require more information, please call our sales team on Tel: +44 (0) 1924 614567













Apollo Healthcare Technologies Limited, Unit A1, Gildersome Spur, Leeds, West Yorkshire, LS27 7JZ Tel: +44 (0) 1924 614567 Fax: +44 (0) 1924 607480 Email: sales@apollo-ht.co.uk www.apollo-ht.co.uk



Saturn Rotate Low[™]

AHSM-078 Rev 1.0 (08.11.24)

The Saturn Rotate Low is an electrically powered rotating profiling chair bed offering independence and mobility to weight bearing users, assisting them through the process of lying down to sitting up and then standing.

The **Saturn Rotate Low** is different to other beds in the Market Place as offers Innovative lateral support, assisting correct posture and aiding fatigue management. It is supplied complete with a mattress designed for users at 'High Risk' of developing pressure ulcers and is an ideal choice for reducing care costs whilst offering independence, dignity and quality of life to the user.















Unique lateral support backrest



High quality motor with low noise level from TiMOTION

Saturn Rotate Low[™]

Featuring all the usual functions of a profiling bed, the **Saturn Rotate Low** also enables the user to rotate and then tilt the chair to aid in the transfer between standing and sitting. The auto-regression function ensures minimal pressure is placed on the sacral region when the bed is profiled, and the unique backrest provides greater lateral support. The Saturn Rotate Low includes a mattress designed for users at 'High Risk' of developing pressure ulcers. It features a vapour permeable, multi-stretch, high frequency welded cover, offering outstanding pressure reduction, comfort and ease of cleaning. A White Glove Delivery and Installation Service is available upon request.

Innovative lateral support backrest assists the upper body to help maintain a safe midline position and also assists with fatigue management

Electrically operated backrest, knee break, height, rotate and tilt



Auto-regression backrest

Profiled bed rotates through 90° to become chair (both left and right sided rotation is available)

Dual caster with central lock (caster diameter: 7.5 cm)

Optional aluminium side rails with quick release rail system (meet clause 201.9.8.3.3.3 of IEC 60601-2-52)

**should not be used when the bed is rotating

12 buttons hand controls with LED illumination for visually impaired




like a rise and recliner chair, assisting with the transfer between sitting and standing



Includes a mattress designed for users at 'High Risk' of developing pressure ulcers, featuring a vapor permeable, multi-stretch, and high-frequency welded cover for enhanced durability and comfort

Comes complete with a set of bed sheets



Handrails to assist with transfer on and off the bed



Scan to watch the product video